

JumpStart[®] Cricket[®]




CASCADE[®]
The daFO People[®]

Sub-malleolar, UCBL trimline

Key Concept — Provides flexible support for young patients with correctable mild to moderate pronation or supination.

The JumpStart Cricket was developed especially for young children with mild to moderate pronation and supination. Stability and alignment improvement come from a deep heel cup, well-defined plantar surface, and medial / lateral forefoot containment in combination with a stable outer foundation.

It requires no casting; the Cricket is fit from size measurements. The Cricket's UCBL trimline allows the orthosis to become enclosed entirely within a shoe. The patient can then practice walking and standing around a centrally-aligned and supported foot position to experience improved stability, balance, and proprioceptive input.



Key Features

Partial-wrap inner liner

cups heel, supports arches, and contains midfoot. Trimline is below top edge of shoe. Made from soft, thin polyethylene plastic.

Solid outer foundation

provides additional support and heel stabilization. Two material options: polyethylene for moderate flexibility or co-poly for firm support.

Prefabricated, fit to measure

Widths: Wide and Narrow

Sizes: 4.00 – 9.00 in. (0.25 in. increments)

Indications

- Fully-correctable foot position
- Mild to moderate low tone pronation or supination

Contraindications

- Uncorrectable foot positions
- Strong pronation or supination

Support Provided

Hindfoot

Full heel cup allows maximum control possible in a low, show insert, trimline.

Forefoot

Partial wrap enclosure allows maximum control possible in a low trimline, helping to maintain heel and forefoot alignment.

Ankle

The Cricket provides no direct ankle control.



No casting needed!



Cascade Dafo, Inc.

1360 Sunset Ave
Ferndale, WA 98248

ph: 800.848.7332
fax: 877.856.2160

intl: + 00 1 360 543 9306
www.cascadedrafo.com